

HI 901 • HI 902

# Multiparameter Food and Dairy Titration Systems

Designed to Measure:

**pH • Acidity • Lactic Acid • Chloride • NaCl • SO<sub>2</sub> • Vitamin C • Alkalinity  
Potassium • Calcium • and more**



**DAIRY**  
line

**ISO 9001:2000**  
CERTIFIED

**HANNA**<sup>®</sup>  
instruments  
With Great Products, Come Great Results™  
[www.hannainst.com](http://www.hannainst.com)

HI 901 • HI 902

# Multiparameter Food and Dairy Titration Systems



HI 901 • HI 902

# Multiparameter Food and Dairy Titration Systems

## List of Methods

- A pH value**  
Calibration and handling of electrodes, measuring pH
- B Preparation and titer determination of the most important titrants**  
NaOH, KOH, HCl, H<sub>2</sub>SO<sub>4</sub>, iodine, thiosulfate, permanganate, AgNO<sub>3</sub>, EDTA
- C Drinking water and mineral waters**  
pH value, acid capacity, water hardness, Ca, Mg, chloride, sulfate, sulfide, chlorine, CO<sub>2</sub>, O<sub>2</sub>, oxidizability
- D Milk, Cheese & Dairy products**  
pH value, acidity, chloride, Ca, vitamin C, Kjeldahl N
- E Edible fats and oils**  
Acid number; free fatty acids, hydroxyl value, iodine value, peroxide values, saponification value
- F Cereals, dry pasta**  
pH value, acidity, NaCl, Kjeldahl N/total protein, Ca, Mg
- G Honey, sugar and sweets**  
pH value, free acid, formol number, reducing sugars
- H Soft drinks, lemonades**  
Citric acid/citrate, phosphoric acid, K, total P
- J Fruit and vegetable juices, jams**  
pH value, total acid, vitamin C, SO<sub>2</sub>, chloride, total P, sulfate, Ca, Mg, K, ash alkalinity, Kjeldahl N, formol number, reducing sugars
- K Beer, vinegar, liquor and wine**  
pH value, total acid, CO<sub>2</sub>, SO<sub>2</sub>, vitamin C, ash alkalinity, chloride, sulfate, total esters
- L Coffee, cocoa and chocolate**  
pH value, acidity, chloride, reducing sugars, Kjeldahl N, free fatty acids

- M Sweeteners, gelling and thickening agents**  
Methoxy- and ethoxy groups, cyclamate, saccharin
- N Fruit, vegetable and mushroom preserves**  
Oxalic acid, NaCl, SO<sub>2</sub>
- O Table salt, spices, pickling salt, seasoning salt, herbal salt and flavored salt**  
NaCl, total iodine, fluoride, tricalciumphosphate, nitrite
- P Meat products, meat extracts, consommé cubes, jellied meat, seasonings, soups, sauces**  
NaCl, Kjeldahl N/raw protein, SO<sub>2</sub>
- Q Qualification of the titrator, validation of a titration method**

## Keep an accurate record of analyses!

HANNA's 900 Series titration systems are easily incorporated into any existing GLP data management program:

- Users can easily record all necessary GLP information with every sample including sample identification, company and operator name, date, time, electrode ID codes and calibration information.
- Data can also be transferred using an integral floppy disk drive for communication with a PC or even other titration systems. Special memory cards are not required.
- All test results can be transferred directly to a PC.
- Users can print reports of analyses directly from the titrator using a standard parallel printer.
- An external monitor and keyboard can be attached for added versatility.
- RS 485 port for future expansion.

## Program View Store Print

with HANNA PC Connectivity



### PC Connectivity

The HI 900 series titrators can be connected to a computer in conjunction with HANNA software.



### VGA Display Connectivity

The information shown on the titrator display can be viewed on a standard VGA display via a 15-pin socket.



### Printer Connectivity

The HI 900 series titrators accept a variety of parallel printers for printing of data.



### PC Keyboard Connectivity

An external PC compatible keyboard can be connected.


www.hannainst.com

## Clip-Lock™ Exchangeable Burette System

With Clip-Lock™, it only takes a couple seconds to exchange the reagent burettes to perform a different titration.

HI 901 and HI 902 Automatic Titrators Feature:

- 320 x 240 pixel LCD w/backlight.
- Precise dosing system (accuracy under 0.1% of burette volume).
- Support up to 100 titration methods (standard and user defined).
- Clip Lock™ — change burettes quickly with auto burette recognition
- Dynamic/Linear dosing feature.
- Fixed end point potential or pH.
- Equivalence point detection (first derivative and second derivative).
- The results are displayed directly in the selected units.
- Titration graph can be displayed on-screen & saved.
- User customized reports can be printed, saved on floppy disk or transferred to PC via RS232 interface.
- Reminders for titrant age and standardization expiration.
- Self diagnostic features for peripheral devices including pump, valve, burette and stirrer.



### A Complete Analysis

These instruments perform a complete analysis comprising of sample preparation, dispensing of titrant solution, stirring, measuring and waiting times, recognition of the end point and storing the results. All the parameters that a titration requires are grouped into a *method*.

The titrators are already supplied with a pack of standard methods or you can create your own. Using a floppy disk or connecting the titrator to the HI 900 PC application, methods (standard and user) can be upgraded, stored or deleted.

# HI 901 • HI 902 Multiparameter Wine Analysis Titration Systems



SPECIFICATIONS	mV	pH	Temperature
Range	-2000.0 to 2000.0 mV	-2.000 to 20.000 pH	-5.0 to 105.0°C/23 to 221°F
Resolution	0.1 mV	0.1/0.01/0.001 pH	0.1°C/0.1°F
Accuracy	±0.1 mV (@25°C/77°F)	±0.001 pH (@25°C/77°F)	±0.1°C/±0.2°F (excluding probe error)
Burette Sizes		5, 10, and 25 mL	
Burette Resolution		1/40000	
Display Resolution		0.001 mL	
Dosing Accuracy		±0.1% of full burette volume	
Display		Graphic LCD, 320 x 240 pixel LCD	
Languages		English, Italian, Portuguese	
Methods		10 standard and 90 user-defined methods	
Burette Auto-Detection		Burette size is automatically recognized when inserted into the unit	
Programmable Stirrer		Propeller type, 100-2500 RPM, automatically held within 10% of the set value, resolution 100 rpm	
Flow Rate		User-selectable from 0.1 mL/min to 2 burette-volumes/min	
pH/mV Measurement		Titration will also perform direct pH and mV measurements	
Temperature Compensation		pH measurements are automatically temperature compensated	
pH Calibration		Manual or automatic at 1-5 points with 4 buffer sets or custom buffers	
Potentiometric Titrations		Acid-Base (pH or mV-Mode), Redox, Precipitation, Complexometric, Non-Aqueous, Ion-Selective, Argentometric (in mV-mode only)	
HI 901 Titration Methods		Fixed mV or pH end-point detection & first equivalency point detection (with the 1st or 2nd derivatives)	
HI 902 Titration Methods		Fixed mV or pH end-point detection & multiple equivalency point detection (with the 1st or 2nd derivatives); back titration	
Measurement Units		User specified expression of concentration units to suit specific calculation requirements	
Real Time & Stored Graphs		mV-Volume or pH-Volume titration curve, 1st derivative curve or 2nd derivative curve, in pH-mode or mV-mode; pH/mV values versus time-data-logging results	
Data Storage:		Up to 100 complete titration and pH/mV logging complete reports	
Disk Drive:		Built-in 3.5" floppy disk drive allows storage and transfer of configurations, preprogrammed methods, custom methods, titration reports and bitmap graph files	
Peripherals		Connections for VGA display, PC-keyboard, parallel printer, RS 232 interface, interface for future expansion	
GLP Conformity		Instrumentation data storage and printing capabilities	
Operating Environment		10 to 40°C (50 to 104°F), up to 95% RH	
Storage Environment		-20 to 70°C (-4 to 158°F), up to 95% RH	
Dimensions		Width x Depth x Height = 390 x 350 x 380 mm (15.3 x 13.8 x 14.9 in)	
Weight		approx. 10 kg (22 lb.) with one pump and stirrer assembly	

Preprogrammed Methods Include:			
Name	Dosing Type	Endpoint Mode	Calculation
NaOH Titrant Concentration	Dynamic	mV equivalency, 1st derivative	Titrant concentration by weight
Titrateable Acidity with Initial pH Measurement	Dynamic	Endpoint titration	Sample concentration by volume
I <sub>2</sub> Titrant Concentration	Dynamic	mV equivalency, 1st derivative	Titrant concentration by weight
Free & Total SO <sub>2</sub> Titrations	Dynamic	mV equivalency, 1st derivative	Sample concentration by volume

### ORDERING INFORMATION

**HI 901** is supplied with (1) 25 mL burette, (1) burette driver, overhead stirrer, glass body pH electrode, temperature probe, power adapter, electrode holder, and instructions.

**HI 902** is supplied with (1) 25 mL burette, (2) burette drivers, overhead stirrer, glass body pH electrode, temperature probe, power adapter, electrode holder, and instructions.

For more information or for a distributor near you: **888.815.6422** • [sales@hannainst.com](mailto:sales@hannainst.com)

Authorized Distributor

