

Products for Beer Analysis

RAFT BREWERIES • BREWERIES • BREWPUBS • HOME BREWERS • CRAFT BREWERIES • BREWERIES • BREWPUBS • HOME BREWER

Temperature



Stainless Steel Probe

HI 9043 Thermometer with 20" or 48" Probe

- Thermocouple Thermometer with a Choice of Either **HI 766TR1** (20") or **HI 766TR2** (48") K-Type Thermocouple Probe
- Allows Temperature Readings From Center of Mash Bed
- Greater Accuracy Than Bimetal Gauge Thermometers
- Faster Response Time than a Thermometer in a Thermowell



Brewing Test Basics

Mash pH and Temperature Values

Know optimum temperature and pH ranges for amylase enzymes. Help maximize yield and product consistency.

Beer Haze

In filtered beer, known values in European Brewery Convention or Formazin Turbidity Units help you to maintain brilliance. If the product is unfiltered, a consistent amount of yeast in suspension is important as well.

Water Quality and Mineral Content

Many beer styles are associated with specific concentrations of ions. Test for brewing water minerals of classic beer styles.

pH

Stainless Steel Probe



HI 99151 Portable pH Meter with Electrode

- Waterproof
- Combination pH/Temp Probe
- Protective Rubber Boot

HI 98127 pH Tester

- Waterproof
- Replaceable Electrode Cartridge
- Economical

Built-in Temp Probe



Beer Haze

HI 94124 Portable Haze Meter

- Small Cuvet Sample for Accurate Readings
- 2 or 3 Point Calibration
- Results in Less Than 30 Seconds

Log and Recall Results

PC Compatible



Your Beer Analysis Resource

Water

HI 3818B Multiparameter Test Kit

- Easy Titration or Comparator Tests
- Includes Chloride, Magnesium, Sulfate and Calcium



ISO 9001:2000 CERTIFIED

(800) 426-6287 • beer@hannainst.com • www.hannainst.com/usa

HANNA[®]
instruments
With Great Products, Come Great Results™

P/L BEERANALYSIS/07_1.0 PRINTED IN USA